

same MSDS as Gold.

PRODUCT SPECIFICATIONS

Synthetic Exam Gloves										Latex Exam Gloves				
PRECISION	SHADOW	NPF	WNF	VITEX	VPF	PHANTOM	PLATINUM	GOLD	SILVER	BRONZE				
Material	Acrylonitrile Butadiene Black Textured	Acrylonitrile Butadiene Blue Textured	Acrylonitrile Butadiene White Textured	Polyvinyl Chloride Cream Smooth	Polyvinyl Chloride Clear Smooth	Rubber Latex Black Textured	Rubber Latex Natural Textured	Rubber Latex Natural Textured	Rubber Latex Natural Smooth	Rubber Latex Natural Textured				
Surface Finish	Textured	Textured	Textured	Smooth	Smooth	Textured	Textured	Textured	Smooth	Textured				
Product Number	SHD932	NPF880	WNF810	VTX990	VPF230	—	PLT550	GLD260	SLV330	BRZ640				
X-Small	PCS772	NPF882	WNF812	VTX992	VPF232	PHM912	PLT552	GLD262	SLV332	BRZ642				
Small	PCS775	NPF885	WNF815	VTX995	VPF235	PHM915	PLT555	GLD265	SLV335	BRZ645				
Medium	PCS776	NPF886	WNF816	VTX996	VPF236	PHM916	PLT556	GLD266	SLV336	BRZ646				
Large	PCS778	NPF888	WNF818	VTX998	VPF238	PHM918	PLT558	GLD268	SLV338	BRZ648				
X-Large														
Dimensions														
Minimum Length:	240 mm	240 mm	240 mm	235 mm	235 mm	240 mm	240 mm	240 mm	240 mm	240 mm				
9.45 inch	9.45 inch	9.45 inch	9.45 inch	9.00 inch	9.00 inch	9.45 inch	9.45 inch	9.45 inch	9.45 inch	9.45 inch				
Minimum Thickness:														
Fingers:	0.06 mm	0.12 mm	0.12 mm	0.06 mm	0.06 mm	0.13 mm	0.13 mm	0.13 mm	0.13 mm	.12 mm				
	2.4 mil	4.8 mil	4.8 mil	2.4 mil	2.4 mil	5.18 mil	5.18 mil	5.18 mil	5.18 mil	4.8 mil				
Palm:	0.06 mm	0.10 mm	0.10 mm	0.06 mm	0.06 mm	0.10 mm	0.10 mm	0.10 mm	0.10 mm	0.10 mm				
	2.4 mil	4.0 mil	4.0 mil	2.4 mil	2.4 mil	4.0 mil	4.0 mil	4.0 mil	4.0 mil	4.0 mil				
Physical Properties														
Tensile Strength (MPa)														
Before Aging	14	14-18	14-18	9	11	19-24	19-24	19-24	19-24	18-22				
After Aging	14	17-22	16-22	9	11	16-22	16-22	16-22	16-22	14-18				
Ultimate Elongation (%)														
Before Aging	650	550-650	530-600	300	360	760-880	700-820	760-880	760-880	650				
After Aging	650	600-720	530-650	300	360	720-820	600-780	720-820	720-820	500				
Protein Content	N/A	N/A	N/A	N/A	N/A	≤50 µg	≤50 µg	≤50 µg	≤50 µg	≤50 µg				
Latex Protein														
Residual Powder Content														
ASTM	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg				
AQL	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5				
Country of Origin	Malaysia	Malaysia	Malaysia	Taiwan	China	Malaysia	Malaysia	Malaysia	Malaysia	Malaysia				
Price	\$\$\$	\$\$\$\$	\$\$\$\$	\$\$\$	\$\$\$	\$\$\$	\$\$\$	\$\$\$	\$\$\$	\$\$\$				

Packaging: Gloves are packaged 100 pcs/box by weight, and 10 boxes/case. XL size are packaged 90 pcs/box by weight, and 10 boxes/case. Specifications are subject to change without prior notice.