

# PRODUCT SPECIFICATIONS

Same MSDS as NPF glove

		Latex Exam Gloves										
		Synthetic Exam Gloves					Latex Exam Gloves					
		PRECISION	SHADOW	NPF	WNF	VITEK	VPF	PHANTOM	PLATINUM	GOLD	SILVER	BRONZE
<b>Material</b>		Acrylonitrile Butadiene	Acrylonitrile Butadiene	Acrylonitrile Butadiene	Acrylonitrile Butadiene	Polyvinyl Chloride	Polyvinyl Chloride	Rubber Latex	Rubber Latex	Rubber Latex	Rubber Latex	Rubber Latex
<b>Color</b>		Violet	Black	Blue	White	Cream	Clear	Black	Natural	Natural	Natural	Natural
<b>Surface Finish</b>		Textured	Textured	Textured	Textured	Smooth	Smooth	Textured	Textured	Textured	Smooth	Textured
<b>Product Number</b>												
<b>X-Small</b>		PCS770	-	NPF880	WNF810	VTX990	VPF230	-	PLT550	GLD260	SLV330	BRZ640
<b>Small</b>		PCS772	SHD932	NPF882	WNF812	VTX992	VPF232	PHM912	PLT552	GLD262	SLV332	BRZ642
<b>Medium</b>		PCS775	SHD935	NPF885	WNF815	VTX995	VPF235	PHM915	PLT555	GLD265	SLV335	BRZ645
<b>Large</b>		PCS776	SHD936	NPF886	WNF816	VTX996	VPF236	PHM916	PLT556	GLD266	SLV336	BRZ646
<b>X-Large</b>		PCS778	SHD938	NPF888	WNF818	VTX998	VPF238	PHM918	PLT558	GLD268	SLV338	BRZ648
<b>Minimum Length:</b>		240 mm	240 mm	240 mm	240 mm	235 mm	235 mm	240 mm	240 mm	240 mm	240 mm	240 mm
		9.45 inch	9.45 inch	9.45 inch	9.45 inch	9.00 inch	9.00 inch	9.45 inch	9.45 inch	9.45 inch	9.45 inch	9.45 inch
<b>Minimum Fitzkness</b>												
<b>Fingers:</b>		0.06 mm	0.12 mm	0.12 mm	0.12 mm	0.06 mm	0.06 mm	0.13 mm	0.13 mm	0.13 mm	0.13 mm	.12 mm
		2.4 mil	4.8 mil	4.8 mil	4.8 mil	2.4 mil	2.4 mil	5.18 mil	5.18 mil	5.18 mil	5.18 mil	4.8 mil
<b>Palm:</b>		0.06 mm	0.10 mm	0.10 mm	0.10 mm	0.06 mm	0.06 mm	0.10 mm	0.10 mm	0.10 mm	0.10 mm	0.10 mm
		2.4 mil	4.0 mil	4.0 mil	4.0 mil	2.4 mil	2.4 mil	4.0 mil	4.0 mil	4.0 mil	4.0 mil	4.0 mil
<b>Physical Properties</b>												
<b>Tensile Strength (MPa)</b>												
Before Aging		14	16-24	14-18	14-18	9	11	19-24	19-24	19-24	19-24	18-22
After Aging		14	17-22	16-22	16-22	9	11	16-22	16-22	16-22	16-22	14-18
<b>Ultimate Elongation (%)</b>												
Before Aging		650	550-650	530-600	530-600	300	360	760-880	700-820	760-880	760-880	650
After Aging		650	600-720	530-650	530-650	300	360	720-820	600-780	720-820	720-820	500
<b>Protein Content</b>												
Latex Protein		N/A	N/A	N/A	N/A	N/A	N/A	≤50 µg	≤50 µg	≤50 µg	≤50 µg	≤50 µg
<b>Residual Powder Content</b>												
ASTM		≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg	≤2 mg
AQL		G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5	G II 2.5
<b>Country of Origin</b>		Malaysia	Malaysia	Malaysia	Malaysia	Taiwan	China	Malaysia	Malaysia	Malaysia	Malaysia	Malaysia
<b>Price</b>		\$\$\$	\$\$\$\$	\$\$\$	\$\$\$	\$\$	\$\$	\$\$\$	\$\$	\$\$	\$\$	\$\$

Packaging: Gloves are packaged 100 pcs/box by weight, and 10 boxes/case. XL size are packaged 90 pcs/box by weight, and 10 boxes/case. Specifications are subject to change without prior notice.