


Cranberry Contour Plus Powder Free Nitrile Exam Gloves,  
Ultra Soft

Contour Plus™ is a patented Cranberry product. Our product was awarded with 4-star rating by prominent consultants at Dental Advisor. Special features include moisturized inner coating with USP grade lanolin and Vitamin E. Contour Plus™ is latex-free powder-free, odourless and comes with Pro White color. It is also First Touch™ enhanced. Gloves are packed 100 pieces per dispenser box.

## Series

XS: 3225      S: 3226      M: 3227      L: 3228      XL: 3229

## Chemical Product and Company Identification

Product Name: Powder-Free Nitrile Examination Gloves, Non-Sterile, Pro White

## Physical Data

Dimension	Extra Small	Small	Medium	Large	Extra Large
Length (mm) min.	240	240	240	240	240
Width (mm)	70 - 79	80 - 89	90 - 99	100 - 109	110 - 120
Thickness (mm) min.					
- Finger	0.11	0.11	0.11	0.11	0.11
- Palm	0.09	0.09	0.09	0.09	0.09
- Cuff	0.06	0.06	0.06	0.06	0.06

Physical Appearance  
 Rolled beaded cuff.  
 Gloves are Pro White color.  
 Micro-roughened.  
 Ambidextrous.

## Stability and Reactivity

Tensile Properties	Aged	Unaged
Tensile Strength (MPa), min.	14	14
Ultimate Elongation (%), min.	400	500

## Ecological Information

Bio-Compatibility: The chemical formulation of the gloves and surface lubricating materials do not contain any substances normally known to be harmful to the user or to any person with whom the gloves contact.

## Other information

To prolong product life, gloves should be stored in a cool, dry location, located away from heat and light.

The data in this Material Safety Data Sheet relates to the specific material designated herein and does not relate to the use in combination with other material or in any process.