


Cranberry Original Lightly Powdered Latex Exam Gloves



Original™ is lightly powdered with low modulus latex. The gloves have a smooth surface and comes in Natural Latex color. It is also First Touch™ enhanced. Gloves are packed 100 pieces per dispenser box.

Series

XS: 3835 S: 3836 M: 3837 L: 3838 XL: 3839

Chemical Product and Company Identification

Product Name: Lightly Powdered Latex Examination Gloves, Non-Sterile, Natural Latex Color

Physical Data

Dimension	Extra Small	Small	Medium	Large	Extra Large
Length (mm) min.	240	240	240	240	240
Width (mm)	70 - 79	80 - 89	90 - 99	100 - 109	110 - 120
Thickness (mm) min.					
- Finger	0.12	0.12	0.12	0.12	0.12
- Palm	0.10	0.10	0.10	0.10	0.10
- Cuff	0.8	0.8	0.8	0.8	0.8

Physical Appearance
 Rolled beaded cuff.
 Gloves are in natural latex color.
 Smooth.
 Ambidextrous.

Stability and Reactivity

Tensile Properties	Aged	Unaged
Tensile Strength (MPa), min.	14	18
Ultimate Elongation (%), min.	500	650

Ecological Information

Bio-Compatibility: The chemical formulation of the gloves and surface lubricating materials do not contain any substances normally known to be harmful to the user or to any person with whom the gloves contact.

Other information

To prolong product life, gloves should be stored in a cool, dry location, located away from heat and light.

The data in this Material Safety Data Sheet relates to the specific material designated herein and does not relate to the use in combination with other material or in any process.