


Cranberry Silkcare Powder Free Examination Gloves



Silkcare™ is a patented Cranberry product. Our product was awarded with 4-Star rating by prominent consultants at Dental Advisor. Special features include moisturized inner coating with USP grade lanolin and Vitamin E. Silkcare™ is powder-free, odourless and comes in Pro White color. It is also First Touch™ enhanced. Gloves are packed 100 pieces per dispenser box.

Series

XS: 7815 S: 7816 M: 7817 L: 7818 XL: 7819

Chemical Product and Company Identification

Product Name: Powder-Free Latex Examination Gloves, Non-Sterile, Natural Latex color

Physical Data

| Dimension | Extra Small | Small | Medium | Large | Extra Large |
|---------------------|-------------|---------|---------|-----------|-------------|
| Length (mm) min. | 240 | 240 | 240 | 240 | 240 |
| Width (mm) | 70 - 79 | 80 - 89 | 90 - 99 | 100 - 109 | 110 - 120 |
| Thickness (mm) min. | | | | | |
| - Finger | 0.17 | 0.17 | 0.17 | 0.17 | 0.17 |
| - Palm | 0.15 | 0.15 | 0.15 | 0.15 | 0.15 |
| - Cuff | 0.10 | 0.10 | 0.10 | 0.10 | 0.10 |

Physical Appearance
 Rolled beaded cuff.
 Gloves are Natural Latex color.
 Micro-roughened.
 Ambidextrous.

Stability and Reactivity

| Tensile Properties | Aged | Unaged |
|-------------------------------|------|--------|
| Tensile Strength (MPa), min. | 16 | 21 |
| Ultimate Elongation (%), min. | 500 | 700 |

Ecological Information

Bio-Compatibility: The chemical formulation of the gloves and surface lubricating materials do not contain any substances normally known to be harmful to the user or to any person with whom the gloves contact.

Other information

To prolong product life, gloves should be stored in a cool, dry location, located away from heat and light.

The data in this Material Safety Data Sheet relates to the specific material designated herein and does not relate to the use in combination with other material or in any process.