

Cranberry Xlim Powder Free Nitrile Examination Gloves, Ultra Soft



Xlim™ is an exclusive Cranberry product. Xlim™ is specially engineered and formulated to provide "2ND SKIN" comfort to dental professionals. It is latex-free, powder-free, odourless and comes in Neo-Blue color. It is also First Touch™ enhanced. Gloves are packed 100 pieces per dispenser box.

Series

XS: 3125 S: 3126 M: 3127 L: 3128 XL: 3129

Chemical Product and Company Identification

Product Name: Powder-Free Nitrile Examination Gloves, Non-Sterile, Neo-Blue

Physical Data

Dimension	Extra Small	Small	Medium	Large	Extra Large
Length (mm) min.	240	240	240	240	240
Width (mm)	70 - 79	80 - 89	90 - 99	100 - 109	110 - 120
Thickness (mm) min.					
- Finger	0.14	0.14	0.14	0.14	0.14
- Palm	0.12	0.12	0.12	0.12	0.12
- Cuff	0.08	0.08	0.08	0.08	0.08
Physical Appearance	Rolled beaded cuff. Gloves are Neo-Blue color. Micro-roughened. Ambidextrous.				

Stability and Reactivity

Tensile Properties	Aged	Unaged
Tensile Strength (MPa), min.	14	14
Ultimate Elongation (%), min.	400	500

Ecological Information

Bio-Compatibility: The chemical formulation of the gloves and surface lubricating materials do not contain any substances normally known to be harmful to the user or to any person with whom the gloves contact.

Other information

To prolong product life, gloves should be stored in a cool, dry location, located away from heat and light.

The data in this Material Safety Data Sheet relates to the specific material designated herein and does not relate to the use in combination with other material or in any process.