



Dealer Inquiries Welcome



Enter Search...

[Home](#) | [Promotions](#) | [About Mydent](#) | [Contact Us](#) | [Redemption Center](#) | [Excelsior](#)

[Home](#) >> [Product Catalog](#) >> [Impression Materials](#) >> [Bite Registration Material](#)

Bite Registration Material

Mousse-like and Thixotropic, fluffy vinyl polysiloxane material. Used for bite registration, copying transfers, and primary impression with bite trays. Thixotropic material has an excellent durometer of over 80 and flows easily, yet will not run or sag. Our "Hard" Bite material #9010 can be used in conjunction with all "laser" CAD/CAM systems.




Size: 2/pk w/6 Mixing Tips

Unit/Case: 10bxs/cs

Weight: 10 lbs

Product Options:

#BR-9001 Unflavored - Regular Set



[DOWNLOAD IMA](#)

Impression Materials Setting Times

Item#	Viscosity	Color	Work Time	Set Time	Total Set Time
VP-9003	Light Body	Regular Set	Yellow	2:30	1:00 3:30
VP-9008		Fast Set	Yellow	1:00	1:00 2:00
VP-9004	Medium Body	Regular Set	Buff	2:30	1:00 3:30
VP-9009		Fast Set	Buff	1:00	1:00 2:00
VP-9005	Monophase Body	Regular Set	Purple	2:30	1:00 3:30
VP-9010		Fast Set	Purple	1:00	1:00 2:00
VP-9007	Heavy Body	Regular Set	Rose	2:30	1:00 3:30
VP-9011		Fast Set	Rose	1:00	1:00 2:00
VP-8006	Putty	Regular Set	Dark Blue	2:00	2:00 4:00
VP-8008		Fast Set	Dark Blue	1:00	2:00 3:00
BR-9001	Bite Registration	Regular Set - Unflavored	Light Blue	2:00	1:00 3:00
BR-9002		Fast Set - Unflavored	Light Blue	:30	1:00 1:30
BR-9003		Super Fast Set - Unflavored	Light Blue	:30	:30 1:00
BR-9004	Registration	Fast Set - Mini Chocolate Chips	Light Blue	:30	1:00 1:30
BR-9005		Super Fast Set - Mini Chocolate Chips	Light Blue	:30	:30 1:00

[Catalog](#) | [Trade Shows](#) | [Rewards Program](#) | [News & PR](#)
[MSDS](#) | [Tattoo Products](#) | [Videos](#) | [Links](#) | [Site Map](#)

2008-2011 © Copyright DEFEND All Rights Reserved