



MATERIAL SAFETY DATA Vista Dental Products

Porcelain Silane

SECTION 1 - PRODUCT IDENTIFICATION

Company Name: Inter-Med, Inc. /Vista Dental Products
2200 Northwestern Avenue
Racine, WI 53404
Phone: (877) 418-4782
Fax: (262) 636-9760
Prepared: 03-31-04

SECTION II - HAZARDOUS INGREDIENTS OF MIXTURES

(ND = NOT DETERMINED) (NA = NOT APPLICABLE)

CAS#	Material	%	OSHA PEL	ACGIH TLV
	Ethyl Alcohol	97		
	Glycidoxypropyltrimethoxysilane	3		

SECTION III - PHYSICAL DATA

Vapor Pressure, mm Hg: ND
Evaporation Rate (ether=1): NE
Solubility in Organic Solvents: Alcohol, Ketones
Specific Gravity (H₂O=1): nd
Odor: Odorless
Vapor Density (Air=1): ND
% Volatile by Volume: ND
Boiling Point: ND
Appearance: Colorless liquid
Melting Point: ND

SECTION IV - FIRE AND EXPLOSION

Flash point: NA
Flammable Limits: NA
Extinguishing Media: Chemical foam, carbon dioxide or dry chemical.
Special Fire Fighting Procedures: None.
Unusual Fire and Explosion Hazards: None.

SECTION V - REACTIVITY DATA

Stability: Stable (x)
Incompatibility: (Materials to avoid) NA
Hazardous Decomposition Products: NA
Hazardous Polymerization: NA
Conditions to Avoid: Extreme heat, keep closed.



**MATERIAL SAFETY DATA Vista Dental Products
Porcelain Silane**

SECTION VI - HEALTH HAZARDS

OSHA Permissible Exposure Limit: None

ACGIH Threshold Exposure Limit: None

Carcinogen, due to the presence of saccharin, Listed in: NTP: OSHA: IARC Monograph:
(X) Not Listed

Other Exposure Limit Used: None

Chronic, Other: None

Acute and Chronic: Irritating to skin and eyes. Avoid ingestion.

Medical Conditions Generally Aggravated by Exposure: None known

Primary Route(s) of Exposure: Eye contact, inhalation and ingestion.

SECTION VII - EMERGENCY AND FIRST AID PROCEDURES

Skin: Wash off well with water immediately.

Eyes: Wash off well with water immediately.

SECTION VIII - SPILL OR LEAK PROCEDURES

Handling & Storage: Avoid extreme heat and light.

Spill Management: Clean up with inert absorbent material and dispose of in an approved manner.

Waste Disposal Methods: Dispose of safely in accordance with local, state and federal regulations.

SECTION IX - PROTECTION INFORMATION/CONTROL MEASURES

Respiratory: Generally not necessary.

Eye Protection: Goggles recommended.

Gloves: Protective gloves recommended.

Ventilation: Generally not necessary.

Other Clothing and Equipment: Lab Coat

Prior to use, the user shall determine the suitability of use for this product.

The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. The information relates only to the specific material designated and does not relate to its use in combination with other materials or its use as to any particular process. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall Inter-Med Inc. / Vista Dental Products be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if Inter-Med Inc. / Vista Dental Products has been advised of the possibility of such damages.