



## MATERIAL SAFETY DATA Vista Dental Products

### Quick Stat™ FS

#### SECTION 1 - PRODUCT IDENTIFICATION

Company Name: Inter-Med, Inc. /Vista Dental Products  
2200 Northwestern Avenue  
Racine, WI 53404  
Phone: (877) 418-4782  
Fax: (262) 636-9760  
Prepared: June 24, 2003

#### SECTION II - HAZARDOUS INGREDIENTS OF MIXTURES

(ND = NOT DETERMINED) (NA = NOT APPLICABLE)

Material	%	OSHA PEL	ACGIH TLV
Ferric Sulfate U.S.P. LD50 (Orally) Mice: 3000 mg/kg	25.0	ND	ND
Amorphous Silica	2.5	ND	ND

Non-hazardous Ingredients

Total: 100

#### SECTION III - PHYSICAL DATA

Vapor Pressure, mm Hg: ND  
Evaporation Rate (ether=1): ND  
Solubility in H<sub>2</sub>O: Highly water soluble  
Specific Gravity (H<sub>2</sub>O=1): ND  
Odor: NA  
Vapor Density (Air=1): ND  
% Volatile by Volume: ND  
Boiling Point: ND  
Appearance: Reddish-brown gel

#### SECTION IV - FIRE AND EXPLOSION

Flash point: ND  
Flammable Limits: NA  
Extinguishing Media: NA  
Special Fire Fighting Procedures: Solution will not burn under normal conditions  
Unusual Fire and Explosion Hazards: None

#### SECTION V - REACTIVITY DATA

Stability: Stable (x)  
Incompatibility: (Materials to avoid) None known  
Hazardous Decomposition Products: None  
Hazardous Polymerization: Will not occur (x)  
Conditions to Avoid: Extreme Heat



**MATERIAL SAFETY DATA Vista Dental Products  
Quick Stat™ FS**

**SECTION VI - HEALTH HAZARDS**

OSHA Permissible Exposure Limit: None

ACGIH Threshold Exposure Limit: None

Carcinogen, due to the presence of saccharin, Listed in: NTP: OSHA: IARC Monograph:  
(X) Not Listed

Other Exposure Limit Used: None

Chronic, Other: None

Acute and chronic: Irritation to eyes and skin. Ingestion may cause vomiting, diarrhea, hemolysis, liver damage.

Medical Conditions Generally Aggravated by Exposure: None known

Primary Route(s) of Exposure: Skin, ingestion

**SECTION VII - EMERGENCY AND FIRST AID PROCEDURES**

Skin: Flush well with water.

Ingestion: Induce vomiting and seek immediate medical advice.

Eyes: Flush well with water.

Inhalation: NA.

**SECTION VIII - SPILL OR LEAK PROCEDURES**

Spill Management: Wash contaminated surfaces with soap and water.

Waste Disposal Methods: Dispose of safely in accordance with local, state and federal regulations.

**SECTION IX - PROTECTION INFORMATION/CONTROL MEASURES**

Respiratory: None required

Eye Protection: Safety Goggles

Gloves: Optional

Ventilation: None required

Other Clothing and Equipment: Lab coat to protect clothing

Avoid excessive heat.

Prior to use, the user shall determine the suitability of use for this product.

The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. The information relates only to the specific material designated and does not relate to its use in combination with other materials or its use as to any particular process. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall Inter-Med Inc. / Vista Dental Products be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if Inter-Med Inc. / Vista Dental Products has been advised of the possibility of such damages.